



Artikel-Nr.: MSTIR
Garantie(Jahre): 2
Hersteller: Vernier Science Education
Erscheinungsdatum: 2016-02-05

Microstirrer

CHF 17.66 *

Inhalt: 1 Stück

inkl. MwSt. [zzgl. Versandkosten](#)

In der Regel an Lager, Lieferzeit ca. 1-3 Werktage

Produktinformationen "Microstirrer"

When Vernier introduced the new Drop Counter, they decided to include a Vernier Microstirrer to enhance stirring. Since then, they've had many requests to have this as a separate product.

The Microstirrer takes the place of the stirring bar that is typically used with a laboratory magnetic stirrer. The new Vernier Microstirrer easily attaches to the end of a pH Sensor, Conductivity Probe, Dissolved Oxygen Probe or Ion-Selective Electrode. When the magnetic stirrer plate (required, but not included) is turned on, its magnetic field turns the paddle wheel of the Microstirrer, resulting in enhanced stirring near the electrode tip. The Microstirrer works with any magnetic stirrer.

Weiterführende Links zu "Microstirrer"

- - [Fragen zum Artikel?](#)
- - [Weitere Artikel von Vernier Science Education](#)